

Мультинаконечники и наконечники

Технические характеристики

По вопросам продаж и поддержки обращайтесь:

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

эл.почта: dbc@nt-rt.ru || сайт: <https://deltalab.nt-rt.ru/>



Multitips extra

Available in different capacities.

Graduation in different colors according to volume.

Each 25 ml and 50 ml box contains a free adapter.

Suitable for (non exclusive)

Biohit® eLine® Lite y Pro, Brand® HandyStep electronic, Eppendorf® Repeater Plus, Multipette® Pro, y Multipette® 4780, Nichimate® Stepper (Nichiryo®), Gilson® Distriman, StepMate®, Ripette®, etc.

When used with Multipette® plus and Multipette® pro the dispenser display is not operational.

Non sterile product.



code	description	colour	sterile	case quantity	case weight (Kg)	case volume (m3)
317-314.1	multitips 1,0 ml	yellow	non sterile	100	0,37	0,006
317-315.1	multitips 2,5 ml	green	non sterile	100	0,48	0,006
317-316.1	multitips 5 ml	blue	non sterile	100	0,51	0,006

317-317.1	multitips 10 ml	orange	non sterile	100	0,62	0,006
317-318.1	multitips 125 ml	red	non sterile	100	0,46	0,006
317-319.1	multitips 50 ml	grey	non sterile	100	0,51	0,006



Multitips

Different capacities, standard and sterile models.
Red graduations.

For use (but not limited to): Eppendorf® Multipipette 4780, Brand® Handystep and Handystep Electronic, Nichimate® Stepper, Ripette® and Ripette® GenX, Minilab® 100/101, Stepmate®, Easystep®, Gilson® Distriman®, etc.

Sterile models are Rnase, Dnase and pyrogen free.

*25ml and 50ml models require adapters.


Available models 316A for the volume of 25 ml and 316B for the volume of 50 ml (ask the sales department).

Individually wrapped.



code	description	sterile	case quantity	case weight (Kg)	case volume (m3)
316-313.1	multitips 0,5 ml	no sterile	100	0,37	0,009
316-314.1	multitips 1,25 ml	no sterile	100	0,48	0,009
316-315.1	multitips 2,5 ml	no sterile	100	0,53	0,009
316-316.1	multitips 5 ml	no sterile	100	0,69	0,009
316-317.1	multitips 12,5 ml	no sterile	100	0,87	0,009

316-318.1	multitips 25 ml	no sterile	25	0,48	0,009
316-319.1	multitips 50 ml	no sterile	25	0,51	0,009
316-313.2	multitips 0,5 ml	sterile R	100	0,50	0,009
316-314.2	multitips 1,25 ml	sterile R	100	0,61	0,009
316-315.2	multitips 2,5 ml	sterile R	100	0,662	0,009
316-316.2	multitips 5 ml	sterile R	100	0,82	0,009
316-317.2	multitips 12,5 ml	sterile R	100	1,02	0,009
316-318.2*	multitips 25 ml	sterile R	25	0,48	0,009
316-319.2*	multitips 50 ml	sterile R	25	0,35	0,009



Robotic Tips for Qiagen (300µl and 1,100µl)

Manufactured in clean room class ISO 8.

Racks and tips in conductive polypropylene.

Free of DNase, RNase and endotoxins.

High precision tips. The conductivity allows the system to recognize the filling height and guarantees a minimum immersion of the tip in the liquid, to ensure safe pipetting and dispensing.

The rack has 96 tips.

They are supplied in double package.

Compatible with BEP 2000®, Biorobot®, Elysis Quattro®, Etimax®,
Evolis®, Triturus® and others.

code	description	rack colour	case quantity	case weight (Kg)	case volume (m3)
318000	robotic tip 300 µl	golden	5×2	1,00	0,004
318001	robotic tip 1,100 µl	silver	5×2	1,80	0,0076



Robotic and conductive tips for TECAN

Manufactured in clean room class ISO 8.

Racks and tips in conductive polypropylene.

Free of DNase, RNase and endotoxins.

High precision tips.

The conductivity allows the system to recognize the filling height and guarantees a minimum immersion of the tip in the liquid, to ensure safe pipetting and dispensing.

The rack has 96 tips.

They are supplied in double package.

Compatible with Tecan Genesis® and Tecan Evo Freedom®..

code	description	rack colour	case quantity	case weight (Kg)	case volume (m3)
318006	50 µl conductive tip	yellow	5×2	0,80	0,045
318002	200 µl conductive tip	blue	5×2	0,80	0,045
318003	1000 µl conductive tip	red	5×2	0,80	0,045
318102	200 µl transparent tip	blue	5×2	0,80	0,045
318103	1000 µl transparent tip	red	5×2	0,80	0,045





2 – 10 ml macrotip

Macro tip for Eppendorf®, Biohit® and Gilson® electronic pipettes.
For use (but not limited to) Gilson Electronic 10000®, Biohit® electronic.



code	description	steril	case quantity	case weight (Kg)	case volume (m3)
200014	tips in bulk	no	10X100	7,24	0,087

200014R	25 tips rack	no	4X25	2,25	0,015
200014RF	25 tips rack filter	STERILE R	4X25	2,20	0,015



50 – 1,000 ul tip

Universal tip with collar. Blue or natural colour.



Suitable (not exclusive): Eppendorf®, Gilson®, Biohit®, DSG®, Elkay®, Finnpipette®, Genex®, Jencons®, Nichiryo®. Suitable for Menarini®.

Rack made of polypropylene, autoclavable up to 121 °C.

Code 200012 minimum unit sale: 1 bag.



code	description	sterile	colour	case qty.	case weight	case vol.
200012	tips in bulk	no	blue	10X1.000	7,95	0,085
200012B	single tips individually wrapper	STERILE R	blue	175	0,27	0,004
200012R	96 tips rack	no	blue	10X96	2,36	0,010



100 – 1,000 ul tip

Universal tip with collar. Blue or natural colour.



Suitable (not exclusive): Eppendorf®, Gilson®, DSG®, Elkay®, Finnpipette®, Genex®, Jencons®, Nichiryo®. Suitable for Menarini®.

Rack made of polypropylene, autoclavable up to 121 °C.

Codes 200070 and 327-30 minimum unit sale: 1 bag.



code	description	sterile	colour	case qty.	case weight	case vol.
200070	tips in bulk	no	natural	10X1.000	8,00	0,083
200070R	96 tips rack	no	natural	10X96	2,40	0,009
200070NR	96 tips rack	STERILE R	natural	10X96	2,38	0,010
327-30	filter tips in bulk	no	natural	10X1.000	8,50	0,083
327-30R	96 tips rack filter	STERILE R	natural	10X96	1,99	0,010



0.1 – 10 µl tip

Natural colour universal crystal tip, graduated. Gilson® type.


Long shape, perfect to pipette samples in deep vessels.

Suitable (not exclusive): Gilson®, Socorex®, Eppendorf®, Biohit®, Nichiryo®, Finnpipette®, Dsg®, Elkay®, Genex®, Jencons®, HTL®, Handrop®.

Rack made of autoclavable polypropylene up to 121 °C.



code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)
301-01	Tips in bulk	NO	15 x 1.000	3,00	0,0260
301-01R	96 tips rack	NO	10 x 96	1,80	0,0100
301-01F	Filter tips in bulk	NO	15 x 1.000	2.92	0,0230
301-01FR	96 filter tips rack	STERILE R	10 x 96	1.25	0,0120





0.1 – 10 ul tip

Natural colour universal crystal tip, graduated. Gilson® type.

Rack made of autoclavable polypropylene up to 121 °C.

Suitable (not exclusive): Gilson®, Eppendorf®, Biohit®, Dsg®, Elkay®, Finnpipette®, Genex®, Jencons®, Nichiryo®, HTL® and Handrop®. Code 326-19 minimum unit sale= 1 bag



code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)
200024	tips in bulk	no	20×1.000	3,80	0,026
199084RC	96 tips rack	no	10×96	1,92	0,009
326-19	filter tips in bulk	no	20×1.000	3,70	0,026
327-19	96 tips rack filter	STERILE R	10×96	1,26	0,010





0.5 – 20 ul tip

Natural colour universal crystal tip, bevelled. Eppendorf® type.

Suitable (not exclusive): Eppendorf®, Gilson®, Biohit®, Brand®, Dsg®, Elkay®, Finnpipette®, Jencons®, Nichiryo®, Socorex®, HTL® and Handrop® monocanal.

Rack made of autoclavable polypropylene up to 121 °C.



code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)
301-11	tips in bulk	no	25X1.000	4,80	0,050
301-10	96 tips rack	no	10×96	1,82	0,009
301-09	96 tips rack	STERILE R	10×96	1,80	0,009
326-37	filter tips in bulk	no	25X1.000	5,40	0,035
327-37R	96 tips rack filter	STERILE R	10×96	1,25	0,010



2 – 20 ul tip

Natural colour universal tip. Eppendorf® type.


Available without filter (2-200µl): See our code 200072.

Suitable (not exclusive): Eppendorf®, Biohit®, Finnpipette®, Gilson®, Socorex®, Brand®, HTL®.

Rack made of autoclavable polypropylene up to 121 °C.



code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)
200010FR	96 tips rack filter	STERILE R	10×96	2,00	0,010



2 – 200 ul natural



Natural colour universal tip, bevelled.
Gilson type®.

Suitable (not exclusive): Gilson® P20 P100 P200, Eppendorf®, Biohit®, Dsg®, Elkay®,
Finnpipette®, Genex®, Jencons®, Socorex®, Nichiryo®, HTL®, etc.

Rack made of autoclavable polypropylene up to 121 °C; with transparent lid and blue base.



code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)
200078	tips in bulk	no	25X1.000	9,20	0,071
199078RC	96 tips rack	no	10×96	1,96	0,009
199078N	96 tips rack	STERILE R	10×96	1,96	0,009





2 – 200 ul universal


Universal tip, bevelled. Gilson® type. Tips with filter: natural colour. Tips without filter: yellow colour.

Suitable (not exclusive): Gilson®, Eppendorf®, Dsg®, Elkay®, Finnpipette® (except Multichannel), Genex®, Jencons®, Nichiryo®, HTL®, Socorex® serie Akura .

Rack made of autoclavable (up to 121 °C) polypropylene.



code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)
327-34	tips with filter in bulk	no	25×1.000	9,20	0,083
327-34R	96 filtered tips in rack	STERILE R	10×96	1,45	0,010
200016	tips in bulk	no	25×1.000	8,81	0,072
200016B	single tips individually wrapped	STERILE R	200	0,19	0,002
200016R	96 tips rack	no	10×96	1,94	0,009



1 – 200 ul natural

Universal natural colour tip with 91 mm length.


Suitable (not exclusive): Pipetman®, Eppendorf®, Nichiryo®, Costar Corning®, Oxford®, Volac®, Socorex®, Titertek®, SMI®, Rainin®, Finnpipette®, etc.

Rack made of polypropylene autoclavable up to 121 °C.

Transparent lid and blue base, with 204 tips in 12 x 17 distribution.



code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)
200026	tips in bulk	no	5X1.000	3,50	0,050
327-26	204 tips rack	STERILE R	8X204	2,23	0,014




1 – 200 ul viscous liquids

Natural colour universal tip. Designed to work with viscous fluids like glycerol or thick tissue homogenates. Suitable also for fragile cellular structures, minimizing the risk of tearing or rupture.

Suitable (not exclusive): Eppendorf®, Socorex®, Nichiryo®, Finnpipette®, Biohit®, Oxford®, Pipetman®, Rainin®, Costar Corning®, Volac®.



code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)
327-40	tips in bulk	no	10X1.000	4,20	0,020



2 – 30 ul tip

Natural colour universal filter tip, bevelled.

Gilson® type.

Suitable (not exclusive): Gilson® P20 P100 P200, Eppendorf®, Biohit®, Dsg®, Elkay®, Finnpiquette®, Genex®, Jencons®, Socorex®, Nichiryo®, HTL®, etc.

Rack made of autoclavable polypropylene up to 121 °C.



code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)
327-36	filter tips in bulk	no	25X1.000	9,50	0,071
327-36R	96 tips rack filter	STERILE R	10×96	1,59	0,009



2 – 100 ul tip

Natural colour universal filter tip, bevelled.



Gilson® type.


Suitable (not exclusive): Gilson® P100 P200, Eppendorf®, Biohit®, Dsg®, Elkay®, Finnpipette®, Genex®, Jencons®, Nichiryo®, HTL®, etc.

Rack made of autoclavable polypropylene up to 121 °C.



code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)
327-23	filter tips in bulk	no	25X1.000	9,00	0,071
327-24	96 tips rack filter	STERILE R	10×96	1,96	0,010





5 – 200 (100) ul tip

Natural colour universal tip. Eppendorf® type.

Suitable (not exclusive): Eppendorf®, Gilson®, Brand®, Biohit®, Dsg®, Elkay®, Finnpipette®, Genex®, Jencons®, Kartell®, Nichiryo®, Socorex®, HTL® and Handrop®.

Rack made of autoclavable (up to 121 °C) polypropylene.



code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)
200072	tips in bulk (5-200uL)	no	25X1.000	9,10	0,071
199014RC	96 tips rack (5-200uL)	no	10×96	1,93	0,010
199014N	96 tips rack (5-200uL)	STERILE R	10×96	1,97	0,010
327-33	tips in bulk (5-100uL)	no	25X1.000	10,00	0,071
327-33R	96 tips rack filter (5-100uL)	STERILE R	10×96	1,97	0,010



0.1 – 200 μ l tip

Universal natural colour tip with collar for gel.

Suitable (not exclusive): Pipetman®, Eppendorf®, Titertek®, Biohit®, Socorex®, Nichiryo®, Costar Corning®, Oxford®, Rainin®, Finnpiquette®.

Translucent rack made of autoclavable polypropylene with hinged lid.



code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)
200079	tips in bulk	no	10X1.000	9,20	0,050
200079R	96 tips rack	STERILE R	10×96	1,20	0,010
200079FS	96 tips rack filter (0,1-100uL)	STERILE R	10×96	1,80	0,010



5 – 200 ul Eppendorf® type


Eppendorf® universal tip.

Suitable (not exclusive): Eppendorf®, Gilson®, Brand®, Biohit®, Dsg®, Elkay®, Finnpipette®, Genex®, Jencons®, Nichiryo®, Socorex®, HTL® and Handrop®.

Rack made of autoclavable polypropylene up to 121 °C.



code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)
200009	tips in bulk	no	25X1.000	8,75	0,072
200009B	single tips individually wrapper	STERILE R	200	0,19	0,002
200009R	96 tips rack	no	10X96	1,98	0,009



5 – 200 ul Gilson® type

Yellow universal tip, Gilson® type.

Suitable (not exclusive)

Gilson®, Eppendorf®, Biohit®, Dsg®, Elkay®, Finnpipette®, Genex®, Jencons®, Nichiryo®, HTL® and Handrop®.

Minimum unit sale: 1 bag.



code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)
199083	tips in bulk	no	25X1.000	9,00	0,071



5 – 200 µl MLA® type

Natural colour universal tip, bevelled. MLA® type.

Translucent rack made of autoclavable polypropylene up to 121 °C; with hinged lid.

Tips in bulk: minimum unit sale: 1 bag.



code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)
200017	tips in bulk	no	25X1.000	10,50	0,071



5 – 200 ul Oxford type

Natural colour universal tip, bevelled. Oxford® type.


Suitable for the DiaMed-Micro typing System machine.

Translucent rack made of autoclavable polypropylene up to 121 °C; with hinged lid.

Tips in bulk: minimum unit sale: 1 bag.



code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)
200019	tips in bulk	no	12X1.000	6,90	0,071



100 – 1,000 ul tip

Universal tip with collar. Blue or natural colour.

Suitable (not exclusive): Eppendorf®, Gilson®, Biohit®, DSG®, Elkay®, Finnpipette®, Genex®, Jencons®, Nichiryo®. Suitable for Menarini®.

Rack made of polypropylene, autoclavable up to 121 °C.




Rack without filter: transparent lid and blue base.

STANDARD filter: hinged transparent lid and transparent base.

Codes 200070, 200012 and 327-30 minimum unit sale: 1 bag.



code	description	sterile	colour	case qty.	case weight	case vol.
200070	tips in bulk	no	natural	10X1.000	8,00	0,083
200070R	96 tips rack	no	natural	10X96	2,40	0,009
200070NR	96 tips rack	STERILE R	natural	10X96	2,38	0,010
327-30	filter tips in bulk	no	natural	10X1.000	8,50	0,083
327-30R	96 tips rack filter	STERILE R	natural	10X96	1,99	0,010



100 – 1.000 ul Gilson® type

Gilson® Universal tip.



Suitable (not exclusive): Gilson®, Eppendorf®, Biohit®, DSG®, Elkay®, Finnpipette®, Genex®, Jencons®, Nichiryo®, Socorex®, HTL® y Handrop®.

Rack made of polypropylene, autoclavable up to 121 °C.

Tips in bulk: minimum unit sale: 1 bag.



code	description	sterile	colour	case qty.	case weight (Kg)	case vol. (m3)
200080.1	tips in bulk	no	blue	10X1.000	9,96	0,085
200080.1R	96 tips rack	no	blue	10X96	2,62	0,015
200082	tips in bulk	no	natural	10X1.000	10,80	0,083
200082R	96 tips rack	no	natural	10X96	2,66	0,015
200082NR	96 tips rack	STERILE R	natural	10X96	2,62	0,013
200082F	filter tips in bulk	no	natural	10X1.000	11,20	0,083
327-16	96 tips rack filter	STERILE R	natural	10X96	2,19	0,013



100 – 1.000 ul Eppendorf® type

Universal tip with collar. Natural colour, Eppendorf® type.

Suitable (not exclusive): Eppendorf®, Nichiryo®, Biohit®, Brand®, Dsg®, Gilson®, Elkay®, Finnpiquette®, Genex®, Jencons®, HTL® and Handrop®. Suitable for Menarini®.

Translucent rack made of autoclavable polypropylene up to 121 °C; with hinged lid.

Tips in bulk: minimum unit sale: 1 bag.



code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)
200029	tips in bulk	no	10X1.000	8,10	0,083
200029R	96 tips rack	no	10X96	2,40	0,010
200029F	filter tips in bulk	no	10X1.000	8,10	0,083
200029RF	96 tips rack filter	STERILE R	10X96	2,40	0,013



100 – 1.000 ul Oxford® type

Universal tip. Oxford® type.

For use (but not limited to): Finnpiquette®, Labsystem®.

Minimum unit sale: 1 bag.

code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)
200020	tips in bulk	no	8×1.000	10,40	0,071




100 – 1.000 ul graduated

Natural colour universal tip, bevelled and graduated.

For use (but not limited to): Gilson®, Rainin®, Finnpiquette®, Nichiryo®, Eppendorf®, Biometra®.

Minimum unit sale: 1 bag.

code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)
200022	tips in bulk	no	10×1.000	10,90	0,080



2 – 300 ul tip

Natural colour universal tip. Biohit® type.


Suitable (not exclusive): Gilson® P200, Finnpipette® (old model), Labsystems®, Nichiryo®, Handrop®, etc.

Rack made of autoclavable polypropylene up to 121 °C.

Tips in bulk: minimum unit sale: 1 bag.



code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)
200074	tips in bulk	no	24X1.000	8,30	0,067
199076RC	96 tips rack	no	10X96	1,82	0,010



2 – 200 µl tip

Yellow universal tip with crown and with a 50 µl margin to prevent any accidental overaspiration.

Suitable (not exclusive): Eppendorf®, Socorex®, Gilson® except NEO model, Nichiryo®, Biohit®, Brand®, Jencons®, Thermo®, DSG®, Elkay®, HTL®, Finnpipette®, Genex®.

Translucent rack made of autoclavable polypropylene up to 121 °C; with hinged lid.

Tips in bulk: minimum unit sale: 1 bag.



code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)
200076	tips in bulk	no	24X1.000	8,82	0,071
200076R	tips in bulk	no	10X96	1,98	0,010
200076RF	96 tips rack filter	STERILE R	10X96	2,04	0,010



20 – 300 ul tip

Natural colour universal tip. Eppendorf® type.



Suitable (not exclusive): Eppendorf®, Gilson®, Brand®, Biohit®, Elkay®, Jencons®, Finnpipette®, Nichiryo®, Socorex® and Handrop®.

Rack made of autoclavable (up to 121 °C) polypropylene.

Tips in bulk: minimum unit sale: 1 bag.



code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)
200077	tips in bulk	no	20X1.000	9,70	0,071
200077R	96 tips rack	no	10X96	1,94	0,009
200077F	filter tips in bulk	no	20X1.000	9,90	0,071
200077FR	96 tips rack filter	STERILE R	10X96	1,92	0,009



50 – 1250 μ l tip

(50-1000 μ l filter tips).

Natural colour universal tip.

Suitable (not exclusive): Socorex®, Gilson®, Nichiryo®, Brand®, HTL®, Biohit®, Eppendorf® Research, DSG®, Elkay®, Finnpipette®, Handrop®.

Rack made of polypropylene, autoclavable up to 121 °C.

Tips in bulk: minimum unit sale: 1 bag.



code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)
200071	tips in bulk	no	10X1.000	8,50	0,083
200071NR	96 tips rack	no	10X96	2,50	0,010
327-31R	96 tips rack filter	STERILE R	10X96	2,80	0,010



5 – 1,000 ul tip

Universal natural colour tip. Beckmann® type.
For use (but not limited to): Sherwood®, Lancer®.
Minimum unit sale: 1 bag.



code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)
200021	tips in bulk	no	10X1.000	6,90	0,071



200 – 1,000 ul tip

Natural colour universal tip special for MLA®. Bevelled.

Suitable (not exclusive): Finnpiquette® colour.

Minimum unit sale: 1 bag.

code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)
200018	tips in bulk	no	8×1.000	2,75	0,071



1 – 5 ml macrotip Eppendorf® type

Natural colour universal macro tip, with collar.

Eppendorf® type.


For use (but not limited to): Biohit®, Dsg®, Eppendorf®, Genex®, Jencons®, Socorex® Acura 835.

Rack made of polypropylene autoclavable up to 121 °C. Polypropylene semi transparent lid and marble base.

Removable lid.



code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)
200075	tips in bulk	no	8X250	5,90	0,088
199075	50 tips rack	no	4X50	1,98	0,015
200075RF	50 tips rack filter	STERILE R	5X50	2,50	0,015



1 – 5 ml blue macrotip Socorex® type

Blue colour universal macro tip with collar.

Socorex® type.

For use (but not limited to): Gilson®, Nichiryo®.


Rack made of polypropylene autoclavable up to 121 °C.

Semi transparent polypropylene lid and marble base.

Removing lid.



code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)
200083	tips in bulk	no	8X250	5,78	0,085
200083R	50 tips rack	no	4X50	2,10	0,015



1 – 5 ml macrotip Finnipipette® type

Natural colour universal macro tip with collar.

Finnipipette® type.

For use (but not limited to): Brand®, Labsystem®.

Rack made of polypropylene autoclavable up to 121 °C.

Semi transparent polypropylene lid and marble base.

Removable lid.



code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)
200028	tips in bulk	no	8X250	6,76	0,085
200028R	50 tips rack	no	4X50	2,20	0,015
200028RF	50 tips rack filter	STERILE R	4X50	2,50	0,015



2 – 10 ml macrotip

Natural colour universal macro tip with collar. Gilson® type.

For use (but not limited to): Labsystem®, Socorex®, Finpipette®, Nichipet®.

Rack made of polypropylene autoclavable up to 121 °C.

Transparent lid and marble base. Sliding lid.



code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)
200015	tips in bulk	no	10X100	7,24	0,087
200015R	25 tips rack	no	4X25	2,25	0,015
200015RF	25 tips rack filter	STERILE R	4X25	2,20	0,015

По вопросам продаж и поддержки обращайтесь:

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

эл.почта: dbc@nt-rt.ru || сайт: <https://deltalab.nt-rt.ru/>