

Воронки

Технические характеристики

По вопросам продаж и поддержки обращайтесь:

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

эл.почта: dbc@nt-rt.ru || сайт: <https://deltalab.nt-rt.ru/>



Funnel holders

Made of autoclavable polypropylene.

Double holder, for funnels from 25 mm (adaptor supplied) up to 170 mm.

Without adaptor, holders can be used with powder funnels with diameter not exceeding 35 mm.

Suitable for use on Ø 8 to 14 mm rods (which include codes 19266 and 19267 for example) by means of adaptors (included).



code	description	case qty.	case weight (Kg)	case vol. (m3)
19269	for 2 funnels	5	0,65	0,001



Rectangular base for burettes

Base made of autoclavable polypropylene, with chromium plated steel rod.

Rods can be screwed one onto the other to obtain different heights.

All models come with 2 rods of Ø 250 x 12 mm.

2 models available: base with central hole, and base with off centre hole.

code	description	presentation	case qty.	case weight (Kg)	case vol. (m3)
19266	base with central hole and + 2 rods	1	1	0,75	0,001
19267	base with off-centre hole and + 2 rods	2	1	0,80	0,005



Burette clamps

Made of autoclavable polypropylene.

Clamps can fit onto rods Ø 8 to 14 mm.

Stainless steel centre arm includes a spring and firmly holds in place the burette by means of rubber pads.

Burette clamps do not obscure graduations.

Suitable for all burettes. Base is supplied separately, see details below.



code	description	presentation	case qty.	case weight (Kg)	case vol. (m3)
19139	for 1 burette	1	5	0,22	0,006
19140	for 2 burette	2	5	0,32	0,013



Büchner funnels

Material: polypropylene.

Unbreakable and resistant to corrosion and heat.

Easily dismantled for cleaning. Autoclavable.



code	capacity	mouth diameter	lower diameter	higher diameter	case qty.	case weight (Kg)	case vol. (m3)
19438	70 ml	55 mm	10,7 mm	15,9 mm	1	0,05	0,007
19440	285 ml	80 mm	11,4 mm	17,7 mm	1	0,13	0,002
19441	390 ml	90 mm	11,2 mm	17,7 mm	1	0,20	0,003
19442	810 ml	110 mm	22,6 mm	30,1 mm	1	0,30	0,004
19443	2.100 ml	160 mm	19,7 mm	30,4 mm	1	0,60	0,012
19445	6.000 ml	240 mm	23,7 mm	37,7 mm	1	1,00	0,028



Rack for sedimentation cones

Material: PMMA.

For 2 cones. Resistant to temperatures from $-40\text{ }^{\circ}\text{C}$ to $85\text{ }^{\circ}\text{C}$.



code	dimensions	case qty.	case weight (Kg)	case vol. (m3)
191056	150 x 300 x 290 mm	1	1,20	0,001



«IMHOFF» sedimentation cone

Material: SAN. Moulded graduations according to DIN standards 12672.

Resistant to temperatures from -20 °C to 85 °C.

Graduation to 1, 2, 4, 6, 8, 10, 20, 30, 40, 60, 80, 100, 200, 300, 400, 500, 600, 700, 800, 900 et 1000 ml.

code	total capacity	mouth diameter	height	case qty.	case weight (Kg)	case vol. (m3)
191055	1.000 ml	125 mm	480 mm	4	1,75	0,019



Powder funnels

Material: polypropylene. Autoclavable. 60° angle.



code	mouth diameter	neck diameter	case qty.	case weight (Kg)	case vol. (m3)
19167	80 mm	15 mm	20	0,28	0,013
19168	100 mm	25 mm	20	0,57	0,011
19169	120 mm	30 mm	1	0,05	0,001
19166	180 mm	40 mm	1	0,09	0,009

[vc_



Large deep trays

Material: polypropylene.
Autoclavable.



Code	External diameter (mm)	Interior diameter (mm)	exit diameter (mm)	length (mm)	neck length (mm)	box quantity	weight (kg)	vo
194171	80	75	5	83	30	1	0.03	0
194172	100	90	8	100	35	1	0.03	0
194173	120	110	9	115	40	1	0.04	0
194175	180	175	12	170	60	1	0.09	0
194176	220	210	17	205	75	1	0.17	0
194177	260	250	21	245	85	1	0.29	0
194178	300	290	24	285	105	1	0.37	0




Long-stem funnels

Material: polypropylene. 60°C angle. Autoclavable.



code	mouth diameter	neck diameter	length	case qty.	case weight (Kg)	case vol. (m3)
19164	80 mm	8 mm	206 mm	10	0,20	0,017



High speed funnels

Material: TPX. With interior helicoid channels for a quicker filtration. Prevents the filter paper from sticking. Transparent. Autoclavable.



code	mouth diameter	neck diameter	length	case qty.	case weight (Kg)	case vol. (m3)
19175	51 mm	3 mm	150 mm	6	0,12	0,0022
19172	70 mm	3 mm	150 mm	6	0,19	0,0042
19173	100 mm	7 mm	108 mm	4	0,26	0,0063
19174	140 mm	10 mm	132 mm	6	0,96	0,0250



Analytical Funnels

Material: polypropylene. Autoclavable.
60° angle.



code	mouth diameter	neck diameter	length	case qty.	case weight (Kg)	case vol. (m3)
19147	45 mm	5 mm	79 mm	20	0,07	0,0036
19148	65 mm	9 mm	112 mm	20	0,20	0,0110
19150	80 mm	10 mm	132 mm	20	0,34	0,0135
19152	100 mm	11 mm	161,5 mm	1	0,03	0,0007
19153	120 mm	11 mm	184 mm	1	0,04	0,0038

По вопросам продаж и поддержки обращайтесь:

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

эл.почта: dbc@nt-rt.ru || сайт: <https://deltalab.nt-rt.ru/>