

Дополнительное оборудование

Технические характеристики

По вопросам продаж и поддержки обращайтесь:

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

эл.почта: dbc@nt-rt.ru || сайт: <https://deltalab.nt-rt.ru/>



Dental containers

Manufactured in polypropylene and intended for the hygienic storage and preservation of dental prostheses, orthodontic sections and dental splints.

It has flexible and waterproof walls and closure system by means of a front flange.

The open container is stackable, which reduces the space of storage, and stackable once closed, thanks to the projections in the base and the lid.

Contact the commercial department for individual packaging.



code	colour	case quantity	case weight (Kg)	case volume (m3)
D01	white	300	9,4	0,070
D02	green	300	9,4	0,070



Digital timers

Digital timer with hours, days and months, 24 h. programme facility, 60 seconds, alarm, includes clip and small bench stand.

These digital timer count forward and down. Both include a LR44 battery.

Code 900400 includes instructions in English and German an code 900600 includes instructions in Spanish, English and French.

code	dimensions	case qty.	case weight (Kg)	case vol. (m3)
900600	64 x 61 x 14 mm	1	0,08	0,0003
900400	52 x 67 x 12 mm	1	0,06	0,0003



One-hand timer

60 Minute times model "Onehand" non-slip surface.
Can be rewound with one hand only.

code	diameter	height	case qty.	case weight (Kg)	case vol. (m3)
19901	72,5 mm	29,5 mm	1	0,24	0,0003



Watch glasses

Made of autoclavable polypropylene. Moulded base for more stability.
Easy cleaning.

code	diameter	unitary weight	case qty.	case weight (Kg)	case vol. (m3)
19454	80 mm	5,5 g	10	0,06	0,002
19455	100 mm	10 g	10	0,10	0,002



Measuring scoops

Material: high density polyethylene.

Made with materials suitable for alimentary use.



code	capacity	length	case qty.	case weight (Kg)	case vol. (m3)
191065	25 ml	135 mm	12	0,08	0,0014
191066	65 ml	185 mm	12	0,20	0,0011
191067	110 ml	215 mm	12	0,35	0,0036
191068	150 ml	250 mm	12	0,42	0,0040
191070	350 ml	310 mm	6	0,39	0,0050
191071	750 ml	350 mm	6	0,69	0,0190
191069	1250 ml	400 mm	6	0,99	0,0115




Measuring scoops

Material: white polypropylene. Autoclavable.

Made with materials suitable for alimentary use.



code	capacity	length	case qty.	case weight (Kg)	case vol. (m3)
191057	10 ml	100 mm	12	0,06	0,0004
191058	25 ml	135 mm	12	0,10	0,0016
191059	50 ml	160 mm	12	0,16	0,0045
191060	100 ml	200 mm	12	0,29	0,0054
191062	250 ml	260 mm	6	0,27	0,0076
191063	500 ml	315 mm	6	0,46	0,0130
191064	1.000 ml	400 mm	6	0,73	0,0150



Laboratory spatulas

Material: polystyrene.

code	description	length	case qty.	case weight (Kg)	case vol. (m3)
19593	spatula-spoon	210 mm	10	0,008	0,0005
19596	spatula-spoon	180 mm	10	0,008	0,0005



Drying rack

PVC coated steel wire drying rack for general labware. It embodies 55 pegs, and 5 brackets specially designed to hold erlenmeyer flasks and beakers. Pegs are coated for a plastic protector* for better protection.

Two rings on the rack allow the rack to be mounted on the wall. It includes a tray in order to avoid drops on the surfaces where the rack is placed.

* The plastic protectors may be red or blue indistinctly.

code	dimensions	case qty.	case weight (Kg)	case vol. (m3)
19212	640 x 360 x 140 mm	1	2,12	0,130



Drying rack

Wall mountable high impact polystyrene drying rack for general labware.

This rack consists of a single mould 4 mm plate with 72 peg sockets.

Push-in pegs can be placed in any configuration on the rack, and are easily removable to accomodate odd shaped items.

Peg sockets have closed ends to avoid dripping and eliminate the risk of biological contamination.

Rack suitable for items with diameter of the mouth larger than 15 mm.

A drip channel collects dripping that drains via a drain tube provided. Includes 11 additional peg sockets (Ø 6 mm) for smaller parts.

Comes complete with a mounting kit and fixing template.

code	dimensions	case qty.	case weight (Kg)	case vol. (m3)
19213	450 x 630 x 110 mm	1	4,00	0,045



Pipette and burette rinsing set

Made of polyethylene and polypropylene.

Can hold pipettes and burettes of length up to 60 cm.

Ensures a fast, safe washing.

Ideal for use in labs where caustic and toxic materials are used.

code	description	case qty.	case weight (Kg)	case vol. (m3)
19215	burette washer	1	3,50	0,11000
19217	pipette washer	1	3,35	0,08000
19218	pipette jar	1	1,80	0,04200
19219	pipette basket	4	2,00	0,02600
19221	pippete jar	1	1,50	0,03000



Vacuum pump

Material: polypropylene. This vacuum pump combines excellent suction capacity with limited water consumption, under low pressure (0.5-1 kg/cm²) as well as high pressure (10 kg/cm²). The built-in, non return valve protects the pump from possible water backflow.

Autoclavable.

Can be dismantled for easy cleaning.



code	description	case qty.	case weight (Kg)	case vol. (m3)
191395	water jet vacuum pump	1	0,03	0,00001



Flask weight ring

Heavy vinyl covered ring ideal to improve stability of Erlenmeyer flasks in water bath.
High chemical resistance.

code	description	int. diameter	case qty.	case weight (Kg)	case vol. (m3)
192470	Erlenmeyer capacity 125 – 500 ml	48 mm	1	0,32	0,0003
192471	Erlenmeyer capacity 250 – 1.000 ml	51 mm	1	0,59	0,0010
192472	Erlenmeyer capacity 500 – 2.000 ml	57 mm	1	0,64	0,0004



Rotary pipette stand

This stand which can hold 94 pipettes consists of 2 discs rotating on a central vertical axis for convenient pipette selection. The lower rotating disc has tapered sides and open bottoms to protect tips and allow water run off. Rotating discs are 18 cm diameter.

The pipette stand comes unassembled with instructions for quick and easy assembly.

Made of autoclavable polypropylene.

Can accomodate all our serology pipettes presented on page 39 and 193, except the 25 ml model.



code	diameter	height	case qty.	case weight (Kg)	case vol. (m3)
19254	230 mm	480 mm	1	0,39	0,005



Pipette stand

Made of polypropylene. Consists of a rounded base with a rod, and a circular pipette stand to hold up to 28 pipettes.

Base and pipette stand are supplied separately.

- Code 19262, rounded base with rod, is specially designed to hold the pipette stand.
- The stand code 19265 accomodates 18 pipettes Ø up to 10 mm, and 10 pipettes Ø up to 15 mm.

The lower disc has little holes for liquid draining. The ring included with the stand allows its raising and lowering on the rod.

code	description	case qty.	case weight (Kg)	case vol. (m3)
19262	rounded base	1	0,24	0,002
19265	18 small and 10 large holes stand	1	0,35	0,007



Filter paper

Our range of filter paper reams is made from high quality cellulose fibres to ensure a good resistance to humidity and a high absorption capacity.

Pore diameter: $\pm$ 30-40 microns. For use in the following applications:

- laboratory worktop protection
- simple filtering operations for varied products
- labware sterilisation
- pharmaceutical product preparation

code	description	grammage	thickness	case qty.	case weight (Kg)	case vol. (m3)
120010	50 x 50 cm sheets	73 g/m2	0,170 mm	1 x 100	2,04	0,006



Absorbent paper

This 125 g/m² absorbent paper provides a total protection thanks to its two layers. The top side consists of a layer of cellulose of great liquid absorption power; the reverse side consists of a layer of polyethylene ensuring that it is completely waterproof.

For use in the following applications:

- laboratory trays and tables protection
- pathological anatomy laboratories
- seed germination (between-paper method)
- spillage recovering using its leaktight surface
- particularly suitable for working with materials that are valuable (precious metals) or hazardous (toxic or corrosive materials)

Paper thickness: 0.250 mm.

Absorption: 240 g/m². Klemm absorption: 86 mm.

code	description	case qty.	case weight (Kg)	case vol. (m3)
120003	50 m x 50 cm reel	1	3,01	0,009
120002	50 x 50 cm sheets	1 x 100	2,90	0,008



Plastic tongue depressors

Made of polystyrene.

Colour: cream. Tongue depressor with smooth rounded edges.

Odour and taste-free. Maintains a high tensile strength.

More flexible than wood depressors.

Features a ribbed finger rest for easy handling. In addition to their traditional use, depressors can be used for spreading balms and ointments, for stirring liquids or mixing medicine, etc.

Size: 16 x 2 cm.

code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)	palet cases
441142.1	-	No	6 x 250	8,30	0,014	108
441142.2	Individual flow pack	STERILE EO	4 x 500	12,66	0,070	24



Wood tongue depressors

Tongue depressor with smooth rounded edges. Made from high-quality wood to provide a smooth, splinter-free finish. Odour and taste-free.

Maintains a high tensile strength.

In addition to their traditional use, depressors can be used for spreading balms and ointments, for stirring liquids or mixing medicine, etc.

The sterile models are supplied with "instructions for use".

Size: 15 x 1.8 cm.

code	description	sterile	case qty.	case weight (Kg)	case vol. (m3)	palet cases
441142	-	No	50 x 100	13,78	0,027	40
441242.2	Individual flow pack	STERILE EO	4 x 500	12,66	0,070	24

По вопросам продаж и поддержки обращайтесь:

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

эл.почта: dbc@nt-rt.ru || сайт: <https://deltalab.nt-rt.ru/>