

Трубки и соединители

Технические характеристики

По вопросам продаж и поддержки обращайтесь:

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35


Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

эл.почта: dbc@nt-rt.ru || сайт: <https://deltalab.nt-rt.ru/>



Flow indicator

Material: SAN. Perfectly transparent.

The rotation of the red ball allows an immediate visual indication of glass or liquid flow.

Very useful for connecting tubes of different diameter from 6.5 to 10 mm.

Not for use with corrosive fluids or temperatures over 85 °C.

code	dimensions	case qty.	case weight (Kg)	case vol. (m3)
19538	102 x 14 x 57 mm	5	0,07	0,0002
19539	90 x 15 x 40 mm	5	0,08	0,0004



Non-return valve

Material: polypropylene.

Autoclavable.

Maximum work pressure: 2 bar (2 kg/cm²).

Minimum work pressure: 0,07 bar, at 20 °C.



code	ext. diameter	case qty.	case weight (Kg)	case vol. (m3)
19418	10,1 / 7,8 mm	10	0,04	0,0002




Composed Disconnectors / connectors

Material: polyethylene. Very useful for connecting rubber or plastic tubing of different diameter.

Composed by two parts.

Connectors fit tightly together.

code	ext. diameter	mouth diameter	case qty.	case weight (Kg)	case vol. (m3)
19435	7,6 / 5,3 mm	2,7 mm	20	0,62	0,0010
19535	12,2 / 9,3 mm	5,5 mm	20	1,06	0,0010



Stopcock, straight bored

Ideal for liquids and gases under low pressure or low vacuum. Material: polypropylene and high density polyethylene that gives a smooth turning action requiring no lubricants.
All are straight bored.

code	mouth diameter	case qty.	case weight (Kg)	case vol. (m3)
19374	4,9 mm	10	0,05	0,0005
19371	7,0 mm	10	0,06	0,0006



Laboratory holders

Ideal for holding tubing or thermometers onto the edge of beakers or other vessels.

Autoclavable.


Minimum order quantity:

19971 – 6 units

19972 – 3 units



code	description	case qty.	case weight (Kg)	case vol. (m3)
19971	for tubing of up to 9.5 mm Ø Polycarbonate	6 x 6	0,13	0,0013
19972	for thermometers from 3.2 mm to 11 mm Ø Polypropylene	12 x 3	0,17	0,0018



«Straight» and 3 way connectors

Material: polypropylene. Autoclavable. Made of one piece.



code	ext. diameter	mouth diameter	case qty.	case weight (Kg)	case vol. (m3)
19510	5,5 / 4,5 / 3 mm	2,0 mm	20	0,16	0,00002
19512	11 / 9 / 6,3 mm	4,5 mm	20	0,06	0,00030
19513	13 / 10,5 / 8,3 mm	6,5 mm	20	0,74	0,00100



«Y» and «+» connectors

Material: polypropylene. Autoclavable.



code	characteristics	ext. diameter	mouth diameter	case qty.	case weight (Kg)	case vol. (m3)
19465	Y	4,7 / 5,7 mm	2,7 mm	20	0,12	0,00018
19467	Y	8,9 / 9,9 mm	5,5 mm	20	0,03	0,00002
19468	Y	10,5 / 11,6 mm	7,3 mm	20	0,09	0,00003



«T» connectors

Material: polypropylene. Autoclavable.



code	characteristics	ext. diameter	mouth diameter	case qty.	case weight (Kg)	case vol. (m3)
19460	T	4,5 / 5,4 mm	2,8 mm	20	0,01	0,0002
19461	T	6,9 / 7,6 mm	4,4 mm	20	0,23	0,0002
19462	T	8,7 / 9,5 mm	6,3 mm	20	0,54	0,0010



Tubing clamps

Material: autoclavable acetal.

These clamps will not distort, rust or corrode.

Different positions to control flow from full open to complete closure.



code	description	case qty.	case weight	case vol.
19505	for tubes of 3.2 to 11 mm (15 positions)	12	0,09	0,0005
19504	for tubes of 11 to 19 mm (12 positions)	6	0,11	0,0007



Silicone tubing

Silicone for numerous applications in laboratories, medicine, pharmacy and industry.

Excellent transparency.

Odourless.

Chemical excellence in general. Excellent durability.

Excellent mechanical properties.

Made of food safe materials.

Latex, phthalate and BPA free material.

Material: High Temperature Vulcanised (HTV) Silicone Elastomer (VMQ).

Service temperature: -40° to 200°C. Tips up to 250°C.

Water-repellent and non-stick.

Gas permeable.

Excellent resistance to UV rays, oxygen, ozone.


Good thermal insulation.



code	int. diameter	ext. diameter	case qty.	case weight (Kg)	case vol. (m3)
350051	0,5 mm	1 mm	1 x 100 m	0,08	0,0011
350013	1 mm	3 mm	1 x 15 m	0,12	0,0010
350024	2 mm	4 mm	1 x 15 m	0,23	0,00 ▲

350034	3 mm	4 mm	1 x 15 m	0,22	0,0043
350035	3 mm	5 mm	1 x 15 m	0,26	0,0026
350037	3 mm	7 mm	1 x 15 m	0,51	0,0032
350046	4 mm	6 mm	1 x 15 m	0,56	0,0026
350048	4 mm	8 mm	1 x 15 m	0,80	0,0043
350057	5 mm	7 mm	1 x 15 m	0,38	0,0039
350059	5 mm	9 mm	1 x 15 m	0,83	0,0047
350069	6 mm	9 mm	1 x 15 m	0,61	0,0100
350610	6 mm	10 mm	1 x 15 m	0,83	0,0096
350612	6 mm	12 mm	1 x 15 m	0,22	0,0043
350079	7 mm	9 mm	1 x 15 m	0,63	0,0091
350710	7 mm	10 mm	1 x 15 m	0,70	0,0076
350810	8 mm	10 mm	1 x 15 m	0,67	0,0091
350812	8 mm	12 mm	1 x 15 m	1,11	0,0054
350814	8 mm	14 mm	1 x 15 m	1,87	0,0008
350912	9 mm	12 mm	1 x 15 m	0,82	0,0120
350915	9 mm	15 mm	1 x 15 m	0,22	0,0043
351014	10 mm	14 mm	1 x 15 m	1,43	0,0081
351420	14 mm	20 mm	1 x 15 m	2,72	0,0240
351622	16 mm	22 mm	1 x 15 m	4,15	0,0320
352030	20 mm	30 mm	1 x 15 m	7,00	0,0490





Silicone sheets

50 x 50 cm silicone sheets ideal as a laboratory benchtop protection or for use in autoclaves to protect any material.

For cleaning simply use a wet cloth.

Resistant up to 260 °C.



code	thickness	case qty.	case weight (Kg)	case vol. (m3)
350002	2 mm	1	0,60	0,030
350004	4 mm	1	1,20	0,001

По вопросам продаж и поддержки обращайтесь:

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

эл.почта: dbc@nt-rt.ru || сайт: <https://deltalab.nt-rt.ru/>